




PADECG consists of professional ECG software and a wireless ECG acquisition module. It provides a mobile solution to meet the demands of caregivers in a mobile treatment environment. EDAN's reliable SEMIP algorithm is also provided in it, which allows you to have a fast and confident ECG assessment before seding the ECG report to Smart ECG Viewer for further analysis.





## Touch Screen

Operating PADECG is just as simple as using your smart phone. Convenient gestures such as long pressing, un-pinching, and dragging are all supported to enable users a better experience.

# ECG Sampling

Both pre-sampling and real-time sampling modes are offered to cater to different situation, no matter it is an emergent case, general ward rounds, or home visiting.

## Filter

Armed with the comprehensive filters, you can acquire superb ECG traces with full confidence just within several touches.

# Algorithm

With EDAN's professional ECG algorithm SEMIP, PADECG gives accurate diagnosis results and offers doctors a reliable reference.

# Minimum Tablet Configuration requirement

	PADECG				
Screen size	7*	Battery	4,000 mAh	os	Android 4.0 OS
Resolution	1024*600 (WSVGA)	Chipset	1.2GHz Dual Core	Memory	16GB + 1GB RAM
Camera	With Auto-focus function	Wi-Fi	Yes	Bluetooth	Yes

# EDAN Edan Instruments, Inc.

- 3/F B, Nanshan Medical Equipments Park, Nanhai Rd 1019#, Shekou, Nanshan Shenzhen, 518067 P.R. China
- Tel: +86-755-26898326
  Fax: +86-755-26898330
  www.edan.com.cn
  Email: info@edan.com.cn
- All rights reserved. Features and specifications are subject to change without notice

